

BRAVO LINE

PEDIATRIC









Model No. -VNP-20

Bravo Line helps find veins at the back of the hands transmission of LED red light, and it can assist medical staff to observe subcutaneous veins of hands, and significantly improve the success rate of the injection.

- Adjustable light, wide application range
- Rechargeable lithium battery.
- Small in size, light in weight, easy to carry.
- Medical cold light, No radiation, safe and reliable.
- Sensitive switch, clear vein and power saving.
- □ Accords with human body engineering, more comfortable grip.

Product Size : 190 x 35 x 35 mm (±1mm)

Product Weight : 84g (±1g)