

## BRV Moon 63+63

Adjustable Colour Temperature



### Special Features

Shadowless Light, Memory Function, High quality sterilizable handle With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

**M** *Memory Function*

**E** *Endo Mode*

### Technical Data

|                        |                                |
|------------------------|--------------------------------|
| Intensity              | 90,000+90,000 Lux ± 10%        |
| Size of Light Field    | 12-30cm                        |
| Colour Temperature     | 4000-5000K                     |
| Colour Rendering Index | 93RA                           |
| LED Life               | >50000Hrs.                     |
| Number of LED          | 63 + 63                        |
| Diameter of Light      | 600mm x 2                      |
| Brightness Control     | Capacitive Touch Panel         |
| Power Supply           | 220V/50Hz AC                   |
| Colour of LED          | White & Yellow (Double Colour) |
| Focusing               | Adjustable                     |
| Battery Backup         | Optional                       |