

## **LED SURGICAL LIGHTS**

## **BRV Moon 63**

Adjustable Colour Temperature



## **Special Features**

Shadowless Light, Memory Function, High quality sterilizable handle With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of head. Body is made of high quality fire resistant polymer.

Mobile Light is having optional feature of power back up. we provide integrated UPS with the unit.





## **Technical Data**

Intensity	90,000Lux ± 10%
Size of Light Field	12-30cm
Colour Temperature	4000-5000K
Colour Rendering Index	93RA
LED Life	>50000Hrs.
Number of LED	63
Diameter of Light	500mm
Brightness Control	Capacitive Touch Panal
Power Supply	220V/50Hz AC
Colour of LED	White & Yellow (Double Colour)
Focusing	Adjustable
Battery Backup	Optional

Easily detachable connector based OT Light / Optional

ture for reference only.

www.bravoplus.in

TRYKER