



This is highly versatile delivery bed offering total flexibility for the doctor, midwife and nurse for the mother it eases greatly the truma of labour and delivery

## Special Features

- Adjustable height and orientation of the seat and back rest offer very high flexibility
- Hand-operated Remote
- U-cut seat for maximum accessibility
- Detachable leg rest meets the requirements of both the professional and the patient
- Leg support for comfortable patient positioning
- Easy- to clean Gynaec Tray
- Firm hand grips provide good support for the patient
- Soft and warm cushion base, with hygienic upper layer
- Foot-operated controls provide operational ease
- Lithotomy
- Adjustable integrated IV pole can be positioned on either side of bed

## Technical Data

Total Length	185 ± 5cm (73")
Length without foot section	140 ± 5cm (55")
Width with side rail	79 ± 5cm (31")
Width without side rail	69 ± 5cm (27")
Minimum Height	69 ± 5cm (19")
Maximum Height	112 ± 5cm (44")
Head section adjustment	-5° to + 70°
Delivery Table Weight	120 kg including all Standard Accessories
Weight capacity	150 kg

## Optional Features

- Hand-operated remote
- Battery back-up supporting up to 50 movements
- Anaesthesia screen holder