


BRV EL- 01

Technical Data

Intensity	40,000Lux ± 10%
Size of Light Field	12-30cm
Colour Temperature	4800K
Colour Rendering Index	93RA
LED Life	>50000Hrs.
Number of LED	14
Brightness Control	Touch
Diameter of Light	150mm
Power Supply	220V/50Hz AC
Battery Backup	Optional



BRV EL- 02

Technical Data

Intensity	25,000Lux ± 10%
Size of Light Field	12-30cm
Colour Temperature	4800K
Colour Rendering Index	93RA
LED Life	>50000Hrs.
Number of LED	7
Diameter of Light	150mm
Power Supply	220V/50Hz AC
Battery Backup	Optional



**" ONE OF
THE BEST MEDICAL
EXAMINATION LIGHT "**