Bravo Plus

FETAL DOPPLER



Model BRV - 201

BRAVO PLUS -201 ultrasonic fetal doppler meets the fetus daily, suitable for routine examination at home, clinic and hospital.

- Easy and convenient to operate
- High- fidelity, crystal clear sound
- Earphone and speaker are feasible
- Various Modes of Operation Waveform, Average, Real time, Manual & Freeze
- User Setable Alarms for FHR
- Display with color LCD
- Unique ergonomic design

TECHNICAL SPECIFICATION

• Ultrasound frequency : 2MHz

• Ultrasound intensity : < 10mW/ cm 2

Power supply
DC Ni-Mh rechargeable battery

AC220/110V, 50/60Hz

• Display : 45mm xX 25mm LCD

Alarm/FHR Measuring range : 50~240 bpm

FHR resolution : 1bpm
FHR accuracy : ± 1bpm
Power consumption : <1W

• Dimension : 135mmxX 95mm xX 35mm

• Weight : 500g