



Model No. : BRV-380

- Lightweight for carrying and Easy - To - Use.
- Silicon Rubber shell protection and stable bracket for table usage.
- Support the oximeter probe for adult, neonate and infant.
- Big size 2.8 inch color TFT LCD display for SPO2/PR/Pulse bargraph plethysmogram.
- Visual and sound alarm Function.
- Adjust the parameters in indicator.
- Low Battery Voltage indicator.
- Automatically Switch off within 3 minutes when no signal.
- Inner Flash memory can store testing result up 300 hours.
- USB interface support upload the data to computer and review the history data with software in PC.
- Standard 4XAAA 1.5V Alkaline Battery or rechargeable Li battery is available for power supply.

SPECIFICATIONS :

SpO2

- Patient : adult, Pediatric , neate
- Range : 35% - 100%
- Resolution : 1%
- Accuracy : 2 % (70% - 100%)
unspecified (less than 70%)

PR

- Range : 25bpm ~ 250bpm
- Resolution : 1bpm
- Accuracy : 1 bpm

TEMPERATURE (OPTIONAL)

- Range : 20 ~ 50° C
- Resolution : 0.1° C
- Accuracy : 0. 2° C
- Dimension : 142X 80 x 29 mm
- Weight : About 90g (with a lithium battery)