

## Model No.: BNJ-06

## **PRODUCTS FEATURES:**

Jaundice is one most common symptoms of the neonatal period, and the clinical Report shows that there is 20 - 30% newborns may appear visibility jaundice The pathologic jaundice for newborn may cause bilirubin encephalopathy, thus find out the pathologic jaundice early as possible is very important. The clinical formerly estimate according the concentration of venous serum bilirubin. In this way, not only waste time and vigor, but also cannot be popularized because of the conditionality. What's more, it can find out the pathologic jaundice and doing the preventive treatment in time, may delay the best treatment time.



The per cutaneous Jaundice tester manufactured by our company can measure the venous serum bilirubin per cutaneous via using optical fiber technology, photoelectric technology, electronic technique and data precessing. It only needs to use the probe of the machine touch the forehead skin of the newborn, and then can measuring date directly and correctly.

The Advantage of our machine are can measuring without ache and without injure, with smaller size and lower weight, easy to operate convenient and safe use, easy to carry and easy promotion to the neonatology departments of all level of hospital for serum bilirubin values per cutaneous dynamic monitoring of the new born. It has the born. It has the broad application prospect in the market.

