

# LABOUR/ DELIVERY/ RECOVERY BED

## Model No. BRV-1003

This newly designed bed enables the obstetrician as well as The lying in woman through all stages of parturition, including labor, delivery and recovery stage with utmost comfort and ease.





## LABOUR/ DELIVERY/ RECOVERY BED



### **Features and Benefits**

- Mattress area divided in to 3 sections made up of CRCA Sheet.
- Detachable polymer moulded head and foot bows provide a sooting look.
- Provided with polymer moulded safety side rails.
- Corner **PU** buffers absorb impact and reduce damage in handling.
- Provided with heavy with duty saline with provision for fixing it at all four locations.
- Castors of 125 /150 mm dia. With central locking mechanism for mobility with control
- Back rest, trendelenburg / reverse trendelenburg highlow, position maneuverable by separate actuator operated by patient control.
- Frame size (approx): 2030L x 900 W X 500-750 H (When fully extended)
- Frame size (approx): 2030L X 900 W X500-750 H (when leg section under the middle section)
- Back rest tilt : 0-75 degree
- Trendelenburg: 0-12 degree
- Reverse trendelenburg: 0-12 degree
- The complete metal components are 7 tank pre treated and epoxy power coated.
- Provided with high density three section mattress.
- Accessories: one collection tray. one set of lithotomy rod and crutches, one heavy duty IV pole, one set of hand grip

#### # Optional

(i) Dust Protective Base Cover

