



B-PAX-1

1

Superior Image Quality

- Optimal image for accurate diagnosis

2

Two Built-in Sensors

- Specialized sensors for Pano & Ceph
- Streamlined workflow and prolonged lifespan of sensors

3

User-friendly Software, EzDent

- It is the complete solution for your consultations

4

Modern and Compact Design

- Raise design value of your clinic
- Save installation space in your practice

Cephalometric (One Shot type)

Superior image quality will be delivered using a Si TFT sensors. Three different ceph image sizes reduce unnecessary X-ray dosage and scans the ideal area of cranial anatomy for your diagnosis and treatment planning.

LATERAL

23X25cm (9X10")



OS (One Shot Standard)

30X25cm (12X10")



OP (One Shot Premium)



ONE SHOT CEPH

With One Shot Ceph, you can acquire high quality images with low X-ray dosage. Plus, short scan time markedly reduces motion artifacts.

PA



Carpus



SMV (Submentovertex)

