



**PATIENT TRANSFER SYSTEM- TRY SWIFT**

We would like to introduce ourselves as suppliers of patient transfer system.

This concept of patient transfer is very Popular in the US, UK and Europe India for the last 3 decades where lifting patients for transferring is not recommended because it causes the nurses and attendants to suffer back pain /back injuries. This has been the protocol followed in all the Emergency Departments around the world including operation Theatres, Wards, ICU'S e.t.c.: basically where ever patient transfer is done between the OT Table to Stretcher or Stretcher to Hospital Bed.

Furthermore it is recommended to use **TRY SWIFT** because the Patient Transfer takes place without lifting the patient's body and is especially useful in the case of accident and emergency cases.

**More than 300 clinics using it for their patient comfort**

**Order Now!**

**FEATURES OF TRY SWIFT**

**TRY SWIFT** is made up of a special poly-carbonate plastic which does break. The length of the board is approximately 170 cms and the width is 50 cms, The Board can be bent over 120 and it will not break,

**FEATURES:**

- Total Support For The Patient's Body From Head To Toe.
- Only 1 Nurse/ Attendant Required For A Successful Patient Transfer.
- No Risk of Back injury & Back Pain To The Nursing Staff who aid in Patient Transfer.
- **TRY SWIFT** is lightweight, Radiolucent & Impact Resistant.
- **TRY SWIFT** made up of special Plastic that is Unbreakable.
- Weights Less Than 4 Kgs, And Can Transfer People With Weight or up to 225 Kgs.
- Lifelong Reusable, If Maintained Properly.
- No Chance of Infection To Patient Or Nurse As It Can Be Cleaned With Any Antiseptic Solution Or Mild Disinfectant After Every Procedure If Required.
- Help To Bridge The Gap & Height Difference Between Trolley And Hospital Bed.
- Worldwide Accepted Compared To The Conventional Rollers Which Hurt The Patient And has A Shelf Life.

Size : 170X50 CM.  
Weight : 3.35 Kgs  
Packing : 1 Pc/ Bag  
Colour : Blue