

**C-ARM COMPATIBLE SEMI-ELECTRIC OT TABLE**



Ref- Page -10  
For Positions

**STEEL QUALITY 304 MEDICAL GRADE STAINLESS STEEL**

**BRV - 2021 EL**

**Special Features**

- Electric up & down by cable remote
- Electric Pillars
- Stainless Fitting
- Detachable, head, leg & Pelvic section
- In Built kidney bridge
- Five sectional Radio - Translucent Top.
- Base & Cover made of attractive, east to clean Stainless Steel.
- Sophisticated mechanics provide smooth, step less articulation of the table Top for precise patient positioning.

**Technical Data**

Length of Table	1900 mm
Width of Table Top	520 mm
Minimum Height (without Mattress)	Min. 750 mm
Maximum Height (without Mattress)	Max. 1000 mm
Trendelenburg	25°
Reverse Trendelenburg	25°
Lateral Tilt	±20°
Back Section	±60°
Head Section	±60°
Leg Section	0°-90°
Floor Locking System	Manual