

DURABLE

ROBUST

UNBREAKABLE JAR

Model No. HVJ-20



FOR LIPO

Main Applications

- Emergency Department
- Intensive Care Unit
- OT / Surgery
- Plastic Surgery
- Obstetrics
- Gastroenterology



TECHNICAL FEATURES

Model No. - HVJ-20	MS	SS	FOR LIPO	
Housing	MS Power Coated on Four castors	Stainless Steel on Four castors	Stainless Steel on Four castors	MS/ Stainless Steel on four castors
Max. Vacuum / LPM	-700 mm Hg, 50 Ltrs / min	-700 mm Hg, 50 Ltrs / min	-700 mm Hg, 70 Ltrs / min	-650 mm Hg, 42 Ltrs / min
Pump Type	Double Rotary Vane Pump	Double Rotary Vane Pump	Double Rotary Vane Pump	Double Rotary Vane Pump
Jars	2X2 Ltrs. PC Jar	2X2 Ltrs. PC Jar	2X2 Ltrs. PC Jar	2X2 Ltrs. PC Jar
Filter	Trap Filter with inbuilt Overflow Protection	Trap Filter with inbuilt Overflow Protection	Trap Filter with inbuilt Overflow Protection	Trap Filter with inbuilt Overflow Protection
Tubing	10 mm ID x 2 mtr. (PVC)	10 mm ID x 2 mtr. (PVC)	10 mm ID x 2 mtr. (PVC)	10 mm ID x 2 mtr. (PVC)
Vacuum Gauge	10 mm dia 0-760 mm Hg	10 mm dia 0-760 mm Hg	10 mm dia 0-760 mm Hg	10 mm dia 0-760 mm Hg
Power Supply	220 V AC, 50 / 60 Hz, 380 Watt	220 V AC, 50 / 60 Hz, 380 Watt	220 V AC, 50 / 60 Hz, 380 Watt	220 V AC, 50 / 60 Hz, 130 Watt
Noise Level	50 dB A ± 3	50 dB A ± 3	50 dB A ± 3	50 dB A ± 3
Dimensions	49 x 36 x 70 cms	49 x 36 x 70 cms	49 x 36 x 70 cms	39 x 31 x 62 cms
Net Weight	37 Kg	36 Kg	36 Kg	25 Kg
Overflow Safety	Dual Mechanical	Dual Mechanical	Dual Mechanical	Dual Mechanical
Optional	Foot Switch / Bacterial Filter / 110V	Foot Switch / Bacterial Filter / 110V	Foot Switch / Bacterial Filter / 110V	Foot Switch / Bacterial Filter / 110V/ Oil- Free Piston Pump

Authorized Certifications from Manufacturer CE, IEC 60601 -1-2012

CE, IEC 60601 -1-2012

CE, IEC 60601 -1-2012

CE, IEC 60601 -1-2012